

Princess[®]

Dare to be perfect



princess[®]



Skin rejuvenation – Wrinkles correction – Facial modeling



Princess Filler... the polyvalent treatment

Indications

Moderate to deep wrinkles, lip contour and lip volume

Injection depth

Mid dermis

Concentration and cross-linking

2.3% HA (23 mg/ml) / BDDE (9%)

Needle

27G


Duration

6 – 9 months





Princess Filler



Forehead furrows

Glabellar lines

Tearthrough

Nasolabial folds

Lip volume and contour,
perioral wrinkles

Marionette lines



Princess Volume... creating volume through a 27G needle

Indications

Deep wrinkles and folds, facial contour and volume

Injection depth

Deep dermis

Concentration and cross-linking

2.3% HA (23 mg/ml) / BDDE (14%)

Needle

27G

Duration

8 to 12 months





Princess Volume



Cheeks

Nasolabial folds

Lip volume

Marionette lines

Facial contour



Princess Rich... the power of HA+glycerol

Indications

Fine lines, mesotherapie

Injection depth

Superficial dermis

Concentration and stabilization

1.8% HA (18 mg/ml) / glycerol (20 mg/ml)

Needle

30G


Treatment

3 sessions (3 weeks interval)



princess[®]

Princess Rich



Face

Neck, décolletage, hands



Why choose Princess?

A unique technology: S.M.A.R.T.

- Based on non-animal HA (biofermentative origin)
- Precise **cross-linking and stabilization** \Rightarrow efficacy, results
- **Viscoelastic optimization** \Rightarrow smooth, homogenous
- Unique **purification** \Rightarrow safety

A product from CROMA

- 30 years of experience – over 15 million syringes sold
- expertise in viscoelastics



Why choose Princess?

Safety

- Full control during production
- Terminal sterilization

	EU standards	Princess Filler	Princess Volume	Princess Rich
BDDE	<2 ppm	<2ppm	<2 ppm	NA
Endotoxins	< 0.5 EU/ml	< 0.5 EU/ml	< 0.5 EU/ml	< 0.5 EU/ml
pH	6,8 – 7,6 ¹	7.0	7.1	7.0
Osmolarity (mosm/kg)	200 - 400 ¹	300	291	300
Proteins	< 1000 ppm ¹²³	11 ppm	8.7 ppm	14 ppm



Why choose Princess?

Amongst the safest products


	Juvederm Ultra	Restylane	Teosyal Global Action	Princess Filler
Concentration (mg/ml)	24	20	25	23
BDDE (ppm)	<2	<2	<2	<2
Endotoxins (EU/ml)	< 0.5	n/a	0.5	< 0.5
Proteins (ppm ¹²³)	2.5	<6	<5	8.7-11



Why choose Princess?

High quality products

- easy and **comfortable to inject** (low injection force)
- immediate, natural and **long-lasting results**
- **highly safe**
- **less pain**
- **less complications**



high satisfaction